## Dilution Rack PBS

1000888

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128 Gefäße / units a 9 mL

## Phosphate Buffered Saline (PBS), sterilized

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

### Preparation of a dilution

(e.g. for the analysis of food specimen)

Open the carton and take out a set of 4 cups. Each of the cups contains a volume of 9 ml PBS. For use of single cups they can be separated at the perforation line.

Make two holes into the sealing foil of the cup. We recommend to use an opener particularly developed for sterile perforation of the aluminium sealing of the device (HyServe Product ID. 1000887).

In case the opener is suspected for bacterial contamination, disinfect the device with cotton wrapped in alcohol and subsequently place it over a flame for sterilisation.

1. After preparing the homogenized specimen, drop 1 mL of the specimen with a steril pipette into the well through one of the holes. (ten times dilution stage,  $10^{-1}$ )
2. Put a pipette into the solution through the other hole and draw up 1 mL. Suck and exhaust the specimen twice to three times for best homogenisation.
3. Take out 1 mL of the diluted homogenized specimen and inoculate onto a Compact Dry plate (see below).
4. In order to generate further dilutions, transfer 1 mL specimen to the next well (dilution  $10^{-2}$ ) and so on. Repeat this step until you get to the dilution calculated.

### Volume

- 9 mL/well

### Formula

- in 1 liter
 

Disodium hydrogen phosphate	0,0785 g
Potassium dihydrogen phosphate	0,0127 g
Sodium chlorid	9,0 g
pH 6,6 – 7,6	

### Storage

- Keep at room temperature (1 – 30 °C)

### Shelf life

- Up to two years after manufacturing.

### Additional Equipment

- Opener, HyServe Product ID 1000887
- Compact Dry Plates (see homepage)