

BagTools®

Essential accessories

A complete range of accessories for a safe and efficient laboratory work.

© 2014 Zeller; a Vötsch Group company
 All rights reserved. No part of this publication may be reproduced, stored in a retrieval system, or transmitted, in any form or by any means, electronic, mechanical, photocopying, recording, or by any information storage and retrieval system, without the prior written permission of Zeller.




BagRack®

Bag holder

BagRack® allows a clean storage of your bags. You can let the bags open for easy pipeting or close it with BagClip®.

3 models: 100 ml, 400 ml and 3500 ml



Option: DripTray® 400 security drip tray in case of liquid leak, for hygiene and cleanliness.



BagOpen®

Bag opener/holder

BagOpen® allows an easy opening of the bag without contact with the sample. Ideal tool for laboratory manipulations.

3 models: 100 ml, 400 ml et 3500 ml




Thanks to BagOpen®, pipeting sample is quick and easy.



BagPipet® & BagTips®

Pipet and sterile tips

BagPipet® pipets fixed chosen volumes (0.1 ou 0.9 ou 1 ml). Used with a BagTips®, one-use sterile straws, easy to eject, BagPipet® allows the pipeting of all products, even viscous, with no blocking risks.




BagPipet® is particularly adapted to viscous and food samples pipeting, after blending in a BagFilter® filter bag with a BagMixer® blender.



BagClip® Closing clip

BagClip® close by pinching any kind of bags, allowing a perfectly tight closing thanks to its three points closing system. BagClip® is trademarked.

3 models: 100 ml, 400 ml et 3500 ml



Thanks to its quick and perfectly tight closing, with BagClip® you can store your samples in a BagRack® in a secured way.



PetriPile® Petri dishes storage stack

PetriPile® allow you to manipulate and transport in a secured manner all your Petri dishes. PetriPile is designed for 36 dishes (16 dishes for the 150 mm PetriPile®).

4 models: Ø 55 / 65 / 90 / 150 mm

Stainless steel designed, stackable and autoclavable, PetriPile® is ideal for incubator loading, according to AFNOR NF V 08-002 standard.

BagTools® Brochure 00311-1 - BagTools®, BagRack®, BagClip®, DripTray®, BagOpen®, BagPipet®, BagTips®, BagFilter®, BagMixer®, et PetriPile® are registered trademarks. Information in this brochure are not contractually binding.



BackRack®
Storage rack
for bags



- Price
- Loan
- Documentation

**BackPipet®
& BackTips®**
Pipet and
sterile tips



- Price
- Loan
- Documentation

BagOpen®
Bag opener/
holder



- Price
- Sample
- Documentation

BagClip®
Closing clip



- Price
- Sample
- Documentation

PetriPile®
Petri dish
storage stack



- Price
- Loan
- Documentation

zeller; a V<

@VUk cf`X'Uh @Vcf[Yf] H'1; 'Ug! F YU[Ybn]Yb
A]_fcV]c' c[]Y! <n]]YbY_cbHfc`Y
-bXi glf]YglfUggY%Z*, () `<c\ YbYa gZ5 i glf]U
HY`Z(' ') +* +* +\$) `: U 'Z(' ') +* +* +\$) +
9a Uj. cZQW4 `UVk cf`X'Uh