




60615-0040



60615-0020



60617-0060



60618-0050


60617-0050

zeller GmbH

Labworld.at
Laborgeräte - Glas - Reagenzien
Mikrobiologie - Hygienekontrolle
Industriestr. 1, 6845 Hohenems, Austria
Tel. +43 (0)5576 76705 Fax +43 (0)5576 76705 7
Email: office@labworld.at


Part No.

- 60617-0050 Dissolver ø 100 mm + shaft
- 60617-0000 Dissolver ø 100 mm
- 60617-0010 Dissolver shaft ø 36 mm
- 60618-0000 Dissolver shaft ø 8x350 mm
- 60617-0060 Dissolver ø 36 mm + shaft



Part No.

- 60618-0050 Stirring blade + shaft
- 60619-0000 Stirring blade (L = 90 mm)
- 60619-0010 Stirring blade (L = 150 mm)
- 60618-0000 Stirring shaft (ø 8x350 mm)
- 60615-0040 Stirring blade tool horizontal (ø 4x230 mm)
- 60615-0020 Stirring blade tool vertical (ø 4x250 mm)



8B00950400



8B00940400



8B00900000



8B00910000



8B00921000

Part No.	Blade Stirrer	Length including blade	Diameter of rod	Diameter of blade
----------	---------------	------------------------	-----------------	-------------------

8B00940400		400 mm	8 mm	70 mm
8B00940800		500 mm	10 mm	70 mm
8B00950400		400 mm	8 mm	70 mm
8B00950800		500 mm	10 mm	70 mm

Part No.	Centrifugal Stirrer	Length including blade	Diameter of rod	Diameter of blade
----------	---------------------	------------------------	-----------------	-------------------

8B00921000		400 mm	8 mm	90/15 mm
------------	--	--------	------	----------

Part No.	Blade Stirrer	Length including blade	Diameter of rod	Diameter of blade
----------	---------------	------------------------	-----------------	-------------------

8B00970000		300 mm	6 mm	50 mm
8B00970300		300 mm	8 mm	50 mm
8B00970400		400 mm	8 mm	50 mm
8B00970500		500 mm	8 mm	50 mm

Part No.	Propeller Stirrer	Length including blade	Diameter of rod	Diameter of blade
----------	-------------------	------------------------	-----------------	-------------------

8B00910000		300 mm	6 mm	50 mm
8B00910300		300 mm	8 mm	50 mm
8B00910400		400 mm	8 mm	50 mm
8B00910500		500 mm	8 mm	50 mm
8B00911000		400 mm	8 mm	70 mm
8B00911400		500 mm	10 mm	100 mm
8B00900000		300 mm	6 mm	50 mm
8B00900400		400 mm	8 mm	50 mm
8B00901000		400 mm	8 mm	100 mm