

Art.Nr.	Typ	Kälteleistung (@20°C)	Kälte- mittel	Arbeits- temperatur	Temperatur- fühler	Temperaturkonstanz DIN 12876	Füll- volumen	Förderstrom max. (0 bar Gegendruck)	Geräusch- pegel	Abmessungen (B x H x T)
0025006624	RC 2 lite	400	R290	-10°C bis RT	PT1000	±0.5 K	1 - 3.5 l	18 l/min	51 dB(A)	225x385x430mm
0004171000	RC 2 basic	400	R134a	-20°C bis RT	PT100	±0.1 K	1.4 - 4.0 l	18 l/min	52 dB(A)	220x475x525mm
0004173000	RC 2 control	400	R134a	-30°C bis RT	PT100	±0.05 K	1.4 - 4.0 l	18 l/min	52 dB(A)	220x475x525mm
0025004186	RC 2 GREEN basic	800	R290	-30°C bis RT	PT100	±0.15 K	1.4 - 4.0 l	21 l/min	61 dB(A)	220x475x525mm
0025006638	RC 2 GREEN control	800	R290	-30°C bis RT	PT100	±0.15 K	1.4 - 4.0 l	21 l/min	61 dB(A)	220x475x525mm
0004181000	RC 5 basic	1400	R290	-30°C bis RT	PT100	±0.2 K	5.0 - 7.0 l	31 l/min	61 dB(A)	310x546x490mm
0004183000	RC 5 control	1400	R290	-30°C bis RT	PT100	±0.1 K	5.0 - 7.0 l	31 l/min	61 dB(A)	310x546x490mm

RC 2 lite
0025006624



RC 2 basic
0004171000



RC 2 control
0004173000



RC 2 GREEN basic
0025004186



RC 2 GREEN control
0025006638



RC 5 basic
0004181000



RC 5 control
0004183000

